

# FISH STOP

## Strip now, fertilize later!




Figure: example of FISH STOP performance for Eurasian perch sperm storage at 2.5°C. Total motility measured by CASA.



Our conservation medium gives you more confort and flexibility in artificial reproduction. Overcome breeders asynchrony, spread breeding over multiple days, transport semen, or dilute very concentrated samples easily with FISH STOP.



Free of antibiotics



Ready to use



Multi-species<sup>1</sup>



Several weeks of storage<sup>2</sup>



FISH STOP is an isotonic culture medium for the dilution and short-term conservation of fish semen at 1-4°C.

<sup>1</sup>: inquire for full-list. <sup>2</sup>: species-specific performance

More info on [www.akymed.com](http://www.akymed.com) or [info@akymed.com](mailto:info@akymed.com)

